

5x Cesium Klentaq AC LA PCR Kit Cat #: 270

Amount: 2 x 1.25 ml (250 reactions)

Shipping conditions: Ice pack

Storage conditions: for best performance, store at -20°C

Shelf life: At least 1 year if stored at -20°C and 10 freeze/thaws or at least 3 months if stored at 4°C.

PRODUCT DESCRIPTION:

Our 5x ready-to-use PCR kit contains Cesium Klentaq AC LA, a cold-sensitive double mutant of Klentaq1 DNA polymerase. Due to its suppressed activity at low temperatures this enzyme is designed for hot-start PCR performance. It also features the Long-and-Accurate feature that allows amplification of longer products with higher fidelity and accuracy. This kit can be used for regular, as well as real-time PCR. It contains everything necessary for a PCR reaction to work perfectly, just add your template, primers/probes and water. For real-time reactions you may need to add a fluorescent dye as an alternative to probes. 5X composition is: 5x Cesium Klentaq AC LA DNA Polymerase, 1 mM dNTPs, 250 mM Tris-Cl pH 9.2, 80 mM ammonium sulfate, 0.125% Brij 58, and 17.5 mM magnesium chloride.

TYPICAL PCR PROTOCOL for a 25 µl reaction:

Reagent	Volume	Final Concentration
5x Cesium Klentaq AC LA PCR Kit reagent	5 µl	1x
Left Primer	variable	0.2 µM
Right Primer	variable	0.2 µM
DNA template [†]	Variable	0.5-100ng
Betaine 5M*	6.5 µl (optional)	1.3 M
de-ionized distilled H ₂ O	Adjust final volume to 25 µl	-

[†] DNA amount depends mostly on genome size and target gene copy number.

*Betaine is a general PCR enhancer. It usually improves the yield and specificity of amplification especially for longer targets.

CYCLING CONDITIONS

1. Pre-incubation: 94° for 2 minutes for 1 cycle
2. Denaturing: 94° for 40-60 seconds
3. Annealing: 55°-70° depending on the specific primers (5° less than T_m) for 40-60 seconds
4. Extension: 68° for 2 min / 1kb target
5. Repeat steps 2-4 for 25-40 cycles

REFERENCES

Kermekchiev, M.B., et al. (2003) Cold-sensitive mutants of Taq DNA polymerase provide a hot start for PCR. Nucl Acids Res. 31, 6139-6147.

Please visit us on the web at www.klentaq.com for troubleshooting and detailed protocols.

Notice to Purchaser

DNA Polymerase Technology products may not be resold, modified for resale or used to manufacture products without an agreement with DNA Polymerase Technology, Inc. Cold sensitive mutant DNA Polymerases by DNA Polymerase Technology are licensed under US Patent No. 6,214,557. No license for Cesium Klentaq AC LA to be used in a Polymerase Chain Reaction has been purchased by DNA Polymerase Technology, Inc.