

## 5x Cesium Klentaq-C PCR Kit

**Amount:** 2 x 1.25 ml (250 reactions)

**Shipping conditions:** Ice pack

**Storage conditions:** for best performance, store at -20°C

**Shelf life:** At least 1 year if stored at -20°C and 10 freeze/thaws or at least 3 months if stored at 4°C.

### PRODUCT DESCRIPTION:

Our 5x ready-to-use PCR kit contains Cesium Klentaq-C, a cold-sensitive double mutant of Klentaq1 polymerase. Due to its suppressed activity at low temperatures this enzyme is designed for hot-start PCR performance. This kit can be used for regular, as well as real-time PCR. It contains everything necessary for a PCR reaction to work perfectly, just add your template, primers/probes and water. For real-time reactions you may need to add a fluorescent dye as an alternative to probes. The 5x Cesium Klentaq-C PCR kit is optimized for targets up to 1 kb in length. **For longer targets, please choose the 5x Cesium Klentaq-C-LA PCR Kit.** 5X composition is: 5x Cesium Klentaq-C DNA Polymerase, 1 mM dNTPs, 250 mM Tris-Cl pH 9.2, 80 mM ammonium sulfate, 0.5% Tween 20, and 17.5 mM magnesium chloride.

### TYPICAL PCR PROTOCOL for a 50µl reaction:

Reagent	Volume	Final Concentration
5x Cesium Klentaq-C PCR Kit reagent	10 µl	1x
Left Primer	variable	0.2 µM
Right Primer	variable	0.2 µM
DNA template <sup>†</sup>	Variable	0.5-100ng
Betaine 5M*	13µl (optional)	1.3 M
de-ionized distilled H <sub>2</sub> O	Adjust final volume to 50µl	-

<sup>†</sup> DNA amount depends mostly on genome size and target gene copy number.

\*Betaine is a general PCR enhancer. It usually improves the yield and specificity of amplification especially for longer targets.

### CYCLING CONDITIONS:

1. Pre-incubation: 94° for 2 minutes for 1 cycle
2. Denaturing: 94° for 40-60 seconds
3. Annealing: 55°-70° depending on the specific primers (5° less than T<sub>m</sub>) for 40-60 seconds
4. Extension: 68° for 2 min / 1kb target
5. Repeat steps 2-4 for 25-40 cycles

### REFERENCES:

Kermekchiev, M.B., et al. (2003) Cold-sensitive mutants of Taq DNA polymerase provide a hot start for PCR. Nucl Acids Res. 31, 6139-6147.

**Please visit us on the web at [www.klentaq.com](http://www.klentaq.com) for troubleshooting and detailed protocols.**

### Notice to Purchaser

DNA Polymerase Technology products may not be resold, modified for resale or used to manufacture products without an agreement with DNA Polymerase Technology, Inc. Cold sensitive mutant DNA Polymerases by DNA Polymerase Technology are licensed under US Patent No. 6,214,557. No license for Cesium Klentaq-C to be used in a Polymerase Chain Reaction has been purchased by DNA Polymerase Technology, Inc.