

## Hot Start Klentaq LA

Cat #: HS110



**Amount:** 100  $\mu$ l enzyme (sufficient for 2000 x 25  $\mu$ l reactions up to 1 kb)

**Shipping conditions:** Ambient

**Storage conditions:** -20°C

**Thermostability:** Retains at least 85% activity after 1 hour at 95°C

**Expiration:** On tube label

### PRODUCT DESCRIPTION:

Hot Start Klentaq LA is made with aptamer-based technology, enabling room temperature reaction set-up.

Klentaq1 is a 5'-exonuclease deficient Taq polymerase (an N-terminal deletion of Taq) with improved fidelity and thermostability. The aptamer binds to the polymerase at sub-cycling temperatures, inactivating the enzyme and preventing spurious amplification. The Long-Accurate feature allows for amplification of longer products with higher fidelity and accuracy. LA enzymes are not recommended for use with dUTP.

10x buffer composition is: 500 mM Tris-Cl pH 9.2, 160 mM ammonium sulfate, 0.5% Brij 58, and 35 mM magnesium chloride. We also offer (upon request) 10x buffer at pH 7.9 for better fidelity.

### TYPICAL PCR PROTOCOL for a 25 $\mu$ l reaction:

Reagent	Volume	Final Concentration
10x Klentaq1 Reaction Buffer	2.5 $\mu$ l	1x
DNTP mix (10 mM)	0.5 $\mu$ l	200 $\mu$ M each
Left Primer	variable	200 nM
Right Primer	variable	200 nM
DNA template†	variable	0.1-100 ng
PCR Enhancer Cocktail (recommended)*	12.5 $\mu$ l	1x
Hot Start Klentaq LA**	0.05 – 0.25 $\mu$ l **	
De-ionized distilled H <sub>2</sub> O	Adjust final volume to 25 $\mu$ l	-

† DNA amount depends mostly on genome size and target gene copy number.

\* For optimal performance, we recommend using one of our PCR Enhancer Cocktails (PEC-1, PEC-1GC, PEC-2, or PEC-2-GC) or 1.3 M Betaine, a general PCR enhancer.

\*\* To determine specific optimal enzyme concentration, we strongly recommend an enzyme titration test for each target. A good starting amount of the enzyme per 25  $\mu$ l reaction is 0.05  $\mu$ l. Targets larger than 1 kb may require more enzyme.

### CYCLING CONDITIONS

1. Denaturing: 94° for 2 minutes for 1 cycle
2. Denaturing: 94° for 40-60 seconds
3. Annealing: 50°-68° depending on the specific primers (5° less than T<sub>m</sub>) for 40-60 seconds
4. Extension: 68° for 2 min/kb target
5. Repeat steps 2-4 for 25-40 cycles

**Please visit us on the web at [www.klentaq.com](http://www.klentaq.com) for troubleshooting and detailed protocols.**

### REFERENCES:

Barnes, W.M. (1994) PCR amplification of up to 35 kb DNA with high fidelity and high yield from bacteriophage templates, PNAS 91, 2216-2220.

U.S. Patent No. 5,436,149

U.S. Patent No. 5,436,149