

Lyoph-Ready Klentaq1

Cat #: GF100

Amount: 4000 x 25 µl reactions up to 1 kb (equivalent to 200ul standard Klentaq1). Volume expansion may occur after dialysis

Shipping conditions: Ice Pack

Storage conditions: 4°C for 4 months or for long-term storage, -20°C for 12 months with up to 10 freeze/thaw cycles

Thermostability: Retains at least 85% activity after 1 hour at 95°C

Expiration: On tube label

PRODUCT DESCRIPTION:

A reduced glycerol preparation of Klentaq1 suitable for lyophilization, is a 5'-exonuclease deficient Taq polymerase (an N-terminal deletion of Taq) with improved fidelity and thermostability. 10x buffer composition is: 500 mM Tris-Cl pH 9.2, 160 mM ammonium sulfate, 0.5% Brij 58, and 35 mM magnesium chloride. We also offer (upon request) 10x buffer at pH 7.9 for better fidelity.

TYPICAL PCR PROTOCOL for a 25 ul reaction:

Reagent	Volume	Final Concentration
10x Klentaq1 reaction buffer	2.5 µl	1x
DNTP mix (10 mM)	0.5 ul	200 uM each
Left Primer	variable	200 nM
Right Primer	variable	200 nM
DNA template†	variable	0.1-100 ng
PCR Enhancer Cocktail (recommended)*	12.5 µl	1x
Klentaq1**	0.05 – 0.25 µl **	
De-ionized distilled H ₂ O	Adjust final volume to 25 ul	-

† DNA amount depends mostly on genome size and target gene copy number.

* For optimal performance, we recommend using one of our PCR Enhancer Cocktails (PEC-1, PEC-1GC, PEC-2, or PEC-2-GC) or 1.3 M Betaine, a general PCR enhancer.

** To determine specific optimal enzyme concentration, we strongly recommend an enzyme titration test for each target. A good starting amount of the enzyme per 25 µl reaction is 0.05 µl. Targets larger than 1 kb may require more enzyme or may benefit from the LA (Long-Accurate) version of the polymerase.

Note: Enzyme amount per reaction should be adjusted based on final enzyme concentration.

CYCLING CONDITIONS

1. Denaturing: 94° for 2 minutes for 1 cycle
2. Denaturing: 94° for 40-60 seconds
3. Annealing: 50°-68° depending on the specific primers (5° less than T_m) for 40-60 seconds
4. Extension: 68° for 2 min/kb target
5. Repeat steps 2-4 for 25-40 cycles

Please visit us on the web at www.klentaq.com for troubleshooting and detailed protocols.

REFERENCES:

Barnes, W.M. (1994) PCR amplification of up to 35 kb DNA with high fidelity and high yield from bacteriophage templates, PNAS 91, 2216-2220.

U.S. Patent No. 5,436,149