



10x RT-LAMP Reaction Buffer

Cat#: RB40

Amount: 2 x 1.5 ml tubes

Shipping conditions: Ambient

Storage conditions: -20° C

Expiration: On tube label

PRODUCT DESCRIPTION:

10x RT-LAMP Reaction Buffer is designed for RT-LAMP with OmniTaq 2.

10x buffer composition is: 500 mM Tris-Cl pH 8.55, 80 mM ammonium sulfate, 0.25% Brij 58, and 29 mM magnesium chloride.

Please visit us on the web at www.klentaq.com for troubleshooting and detailed protocols.