



## **10x RT-LAMP Reaction Buffer**

**Cat#: RB40**

**Amount:** 1.5 ml

**Shipping conditions:** Ambient

**Storage conditions:** -20° C

**Shelf life:** At least 1 year -20° C

### **PRODUCT DESCRIPTION:**

10x RT-LAMP Reaction Buffer is designed for RT-LAMP with OmniTaq 2.

10x buffer composition is: 500 mM Tris-Cl pH 8.55, 80 mM ammonium sulfate, 0.25% Brij 58, and 29 mM magnesium chloride.

**Please visit us on the web at [www.klentaq.com](http://www.klentaq.com) for troubleshooting and detailed protocols.**