



10x Klentaq Mutant Reaction Buffer

Cat#: RB10

Amount: 1.5 ml

Shipping conditions: Ambient

Storage conditions: -20° C

Shelf life: At least 1 year -20° C

PRODUCT DESCRIPTION:

10x Klentaq Mutant Reaction Buffer is specially formulated for optimal performance with our Klentaq Mutant enzymes (Omni Klentaq, Omni Klentaq 2, Cesium Klentaq, etc.)

10x buffer composition is: 500 mM Tris pH 9.2, 160 mM ammonium sulfate, 0.25% Brij 58, and 35 mM magnesium chloride.

Please visit us on the web at www.klentaq.com for troubleshooting and detailed protocols.