



10x Klentaq-S Reaction Buffer

Cat#: RB70

Amount: 2 x 1.5 ml tubes

Shipping conditions: Ambient

Storage conditions: -20° C

Expiration: On tube label

PRODUCT DESCRIPTION:

10x Klentaq-S Reaction Buffer is specially formulated for optimal performance with Klentaq-S DNA polymerase.

10x buffer composition is: 500 mM Tris-Cl pH 7.8, 160 mM ammonium sulfate, 0.25% Brij 58, and 35 mM magnesium chloride. The enzyme may not perform as well at a higher pH.

Please visit us on the web at www.klentaq.com for troubleshooting and detailed protocols.