



10x Klentaq1 Reaction Buffer pH 7.9

Cat#: RB21

Amount: 1.5 ml

Shipping conditions: Ambient

Storage conditions: -20° C

Shelf life: At least 1 year -20° C

PRODUCT DESCRIPTION:

10x Klentaq1 Reaction Buffer pH 7.9 is high-fidelity PCR reaction buffer for use with Klentaq1 and Klentaq LA.

10x buffer composition is: 500 mM Tris-Cl pH 7.9, 160 mM ammonium sulfate, 0.5% Brij 58, and 35 mM magnesium chloride.

Please visit us on the web at www.klentaq.com for troubleshooting and detailed protocols.

Notice to Purchaser DNA Polymerase Technology products may not be resold, modified for resale or used to manufacture products without an agreement with DNA Polymerase Technology, Inc