



10x Klentaq1 Reaction Buffer pH 7.9

Cat#: RB21

Amount: 1.5 ml

Shipping conditions: Ambient

Storage conditions: -20° C

Shelf life: At least 1 year -20° C

PRODUCT DESCRIPTION:

10x Klentaq1 Reaction Buffer pH 7.9 is high-fidelity PCR reaction buffer for use with Klentaq1 and Klentaq LA.

10x buffer composition is: 500 mM Tris-Cl pH 7.9, 160 mM ammonium sulfate, 0.5% Brij 58, and 35 mM magnesium chloride.

Please visit us on the web at www.klentaq.com for troubleshooting and detailed protocols.