



## **10x Klentaq1 Reaction Buffer pH 9.2**

**Cat#: RB20**

**Amount:** 1.5 ml

**Shipping conditions:** Ambient

**Storage conditions:** -20° C

**Shelf life:** At least 1 year -20° C

### **PRODUCT DESCRIPTION:**

10x Klentaq1 Reaction Buffer pH 9.2 is our standard PCR reaction buffer for use with Klentaq1 and Klentaq LA.

10x buffer composition is: 500 mM Tris pH 9.2, 160 mM ammonium sulfate, 0.5% Brij 58, and 35 mM magnesium chloride.

**Please visit us on the web at [www.klentaq.com](http://www.klentaq.com) for troubleshooting and detailed protocols.**