



10x Taq Mutant Reaction Buffer

Cat#: RB30

Amount: 1.5 ml

Shipping conditions: Ambient

Storage conditions: -20° C

Shelf life: At least 1 year -20° C

PRODUCT DESCRIPTION:

10x Taq Mutant Reaction Buffer is specially formulated for optimal performance with our Taq Mutant enzymes (OmniTaq, OmniTaq 2, OmniTaq 3, Cesium Taq, etc.)

10x buffer composition is: 500 mM Tris pH 9.1, 160 mM ammonium sulfate, 0.25% Brij 58, and 25 mM magnesium chloride.

Please visit us on the web at www.klentaq.com for troubleshooting and detailed protocols.

Notice to Purchaser DNA Polymerase Technology products may not be resold, modified for resale or used to manufacture products without an agreement with DNA Polymerase Technology, Inc