



## **10x Taq Mutant Reaction Buffer**

**Cat#: RB30**

**Amount:** 1.5 ml

**Shipping conditions:** Ambient

**Storage conditions:** -20° C

**Shelf life:** At least 1 year -20° C

### **PRODUCT DESCRIPTION:**

10x Taq Mutant Reaction Buffer is specially formulated for optimal performance with our Taq Mutant enzymes (OmniTaq, OmniTaq 2, OmniTaq 3, Cesium Taq, etc.)

10x buffer composition is: 500 mM Tris pH 9.1, 160 mM ammonium sulfate, 0.25% Brij 58, and 25 mM magnesium chloride.

**Please visit us on the web at [www.klentaq.com](http://www.klentaq.com) for troubleshooting and detailed protocols.**