



10x Taq Mutant Reaction Buffer

Cat#: RB30

Amount: 2 x 1.5 ml tubes

Shipping conditions: Ambient

Storage conditions: -20° C

Expiration: On tube label

PRODUCT DESCRIPTION:

10x Taq Mutant Reaction Buffer is specially formulated for optimal performance with our Taq Mutant enzymes (OmniTaq, OmniTaq 2, OmniTaq 3, Cesium Taq, etc.)

10x buffer composition is: 500 mM Tris pH 9.1, 160 mM ammonium sulfate, 0.25% Brij 58, and 25 mM magnesium chloride.

Please visit us on the web at www.klentaq.com for troubleshooting and detailed protocols.